(English)

Universal F®

## Anesthesia Breathing Circuit


Keep these instructions for future reference.

### Intended Use/Indications:

This breathing circuit is a device that is intended to administer medical gases and/or anesthetic gases to a patient during anesthesia for inhalation, or medical gases to a patient during respiratory care inhalation.

### Directions for Use:

Illustration Pertains to Procedure Only.

1. Install inspiratory limb connector (A), containing the blue-green colored inner tube, to the absorbers inspiratory port.
2. Install expiratory connector hose (B) to absorber's expiratory port.
3. **PRIOR TO PATIENT CONNECTION:** Examine circuit for any damage which may have occurred during shipping. Also, inspect the inner tube for kinks. Pressurize the circuit to confirm there is no leakage and/or occlusion of the breathing circuit.
4. Attach the face mask or endotracheal tube to the mask elbow.
5. Discard circuit after single patient use.

**CAUTION:** If the individually packaged product indicates the presence of natural rubber latex, the following statement applies: This Product Contains Natural Rubber Latex Which May Cause Allergic Reactions. Protect from ultraviolet light. Failure to do so may cause deterioration of the latex components.

**CAUTION:** Federal law (U.S.) restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

Test Data available upon request.

Manufacturer

**KINGSYSTEMS CORPORATION**

Noblesville, IN USA 46060 TEL 317-776-6823 FAX 317-776-6827

**EU Regulatory Representative**

Medical Product Service GmbH  
Borngasse 20  
35619 Braunfels, Germany

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