(English)

Universal F2®

Anesthesia Breathing Circuit


Keep these instructions for future reference.

Intended Use/Indications:

This breathing circuit is a device that is intended to administer medical gases and/or anesthetic gases to a patient during anesthesia for inhalation, or medical gases to a patient during respiratory care inhalation.

Directions for Use:

Illustration Pertains to Procedure Only.

1. Install breathing bag on absorber bag mount.
2. Install co-axial limb connector (A), to the King Systems F2® manifold on the anesthesia machine.
3. **PRIOR TO PATIENT CONNECTION:** Examine circuit for any damage which may have occurred during shipping. Also, inspect the inner tube for kinks. Pressurize the circuit to confirm there is no leakage and/or occlusion of the breathing circuit.
4. Attach the face mask or endotracheal tube to the mask elbow.
5. Discard circuit after single patient use.

CAUTION: If the individually packaged product indicates the presence of natural rubber latex, the following statement applies: This Product Contains Natural Rubber Latex Which May Cause Allergic Reactions. Protect from ultraviolet light. Failure to do so may cause deterioration of the latex components.

CAUTION: Federal law (U.S.) restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

Test Data available upon request.

Manufacturer

KINGSYSTEMS CORPORATION

Noblesville, IN USA 46060 TEL 317-776-6823 FAX 317-776-6827

EU Regulatory Representative

Medical Product Service GmbH
Borngasse 20
35619 Braunfels, Germany

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