

# Ambu® VivaSight™ Endobronchial Blocker Tube

## DATASHEET

VivaSight Endobronchial Blocker Tube is a single-use endobronchial blocker for One Lung Ventilation (OLV).

It consists of a “steerable” balloon-tipped catheter that is visually guided to a selected bronchial segment to affect balloon blockage of the distal airway.

The angled distal tip of VivaSight Endobronchial Blocker Tube can be adjusted to facilitate navigation to the bronchi.

### One Lung Ventilation made simple

For visual guidance and positioning the VivaSight Endobronchial Blocker Tube can be used with VivaSight 2 SLT single lumen tube with integrated camera.

When used in conjunction with VivaSight 2 SLT, continuous visualization of the tube and endobronchial blocker throughout the procedure is provided.

For One Lung Ventilation the Ambu visualization portfolio also includes Ambu® VivaSight™ 2 DLT, a double lumen tube with integrated camera, VivaSight 2 SLT, a single lumen tube with integrated camera, as well as the single-use bronchoscope aScope™ 4 Broncho Slim.

When using VivaSight Endobronchial Blocker together with VivaSight 2 SLT, the devices can be connected to the Ambu® aView™ 2 Advance Monitor, Ambu’s full HD monitor for high quality imaging and Dual View.

Alternatively, VivaSight 2 SLT can also be connected with the Ambu® aView™ Monitor.



## KEY BENEFITS

Compatible with Ambu® VivaSight™ 2 SLT for visualization of the SLT and endobronchial blocker during OLV procedure

# Ambu

## SPECIFICATIONS

Ambu® VivaSight™ Endobronchial Blocker Tube	Size	Length (±5 mm)
412900000	9 Fr	695 mm

**General:** VivaSight Endobronchial Blocker Tube features a blue silicone cuff  
**Usage:** Single-use only. Sterile (ETO) unless package has been opened or damaged  
**Pack:** Supplied individually packaged  
**Storage:** Store in a dry, cool and dark place  
**Shelf life:** 5 years  
 This product is not made with phthalates


## ORDERING SPECIFICATIONS

Ambu® VivaSight™ Endobronchial Blocker Tube	Order number	Quantity/box
Ambu® VivaSight™ Endobronchial Blocker Tube	412900000	5

Ambu® VivaSight™ 2 SLT Kits	Order number	Quantity/box
Ambu® VivaSight™ 2 SLT 7.0 mm (with adapter cable)	412701000	5
Ambu® VivaSight™ 2 SLT 7.5 mm (with adapter cable)	412751000	5
Ambu® VivaSight™ 2 SLT 8.0 mm (with adapter cable)	412801000	5

Ambu® VivaSight™ 2 SLT + EBB Kits	Order number	Quantity/box
Ambu® VivaSight™ 2 SLT + EBB Kit - 7.0 mm, 9 Fr (with adapter cable)	412709000	5
Ambu® VivaSight™ 2 SLT + EBB Kit - 7.5 mm, 9 Fr (with adapter cable)	412759000	5
Ambu® VivaSight™ 2 SLT + EBB Kit - 8.0 mm, 9 Fr (with adapter cable)	412809000	5

Monitors	Order number	Quantity/box
Ambu® aView™ 2 Advance Monitor	405011000	1
Ambu® aView™ Monitor	405002000	1



Ambu® VivaSight™ Endobronchial Blocker used with Ambu® VivaSight™ 2 SLT and Ambu® aView™ 2 Advance Monitor

# Ambu

Ambu® VivaSight™  
 Endobronchial Blocker  
 Tube is manufactured by:  
 Well Lead Medical Co., Ltd.  
 C-4 Jinhu Industrial Estate  
 Hualong, 511434 Panyu,  
 Guangzhou  
 People's Republic of China

**EU Importer**  
 Ambu A/S  
 Baltorpbakken 13  
 2750 Ballerup  
 Denmark  
 T +45 7225 2000  
 ambu.com

