

## DATASHEET

CasPeR the CPR Dog incorporates all necessary features for teaching basic dog CPR for pet owners.

It allows training of mouth-to-snout resuscitation and chest compressions with proper hand placement. In addition a manually generated femoral pulse can be checked.

### Sanitary System

The individual mouth/nosepieces are easy to clean and the disposable airway system with a non-re-breathing valve prevents cross-contamination and is changed once a day.

# CasPeR the CPR Dog



## KEY BENEFITS

- Realistic features
- Sanitary system prevents cross-contamination: single-use face shield/lung system Easy and quick change of airways
- Femoral pulse

## ORDERING SPECIFICATIONS

### Item no.

100-5000U	CasPeR the CPR Dog incl. pulse bulb
LF01158U	Disposable Lower Airways, CPR Dog (Pkg. 10)
LF01164U	Disposable Nosepieces, CPR Dog (Pkg. 10)

# Ambu

■ Ambu A/S  
Baltorpbakken 13  
2750 Ballerup  
Denmark  
T +45 72 25 20 00  
F +45 72 25 20 50  
ambu.com