

DATASHEET

Choking emergencies are not selective when it comes to age groups, so several generations of choking manikins are available to choose from.

Each manikin functions on the principle of clearing a blocked airway with forced air generated with the proper procedure; so abdominal thrusts or back blows can be trained.

Therefore, the manikins feature anatomical landmarks including the rib cage, xiphoid process and suprasternal notch.

With the appropriate force, a direct visual confirmation of an expelled object is given to the students.



Choking Manikins



KEY BENEFITS

- 3 generations of choking manikins
- Training of abdominal thrusts or back blows
- Anatomical landmarks
- Visual confirmation of expelled object
- Foreign objects included

ORDERING SPECIFICATIONS

Item no.	
100-1602U	Adult Choking Manikin Adult torso with shirt, 2 foreign objects, soft carry bag
100-1615U	Adolescent Choking Manikin Adolescent torso with shirt, 2 foreign objects, soft carry bag
100-1620U	Child Choking Manikin Child torso with shirt, 2 foreign objects, soft carry bag
100-1640U	Infant Choking Manikin Infant fullbody manikin with pants/shirt, 1 foreign object, soft carry bag
100-1630U	Obese Choker Manikin with Bag Obese mature torso with shirt, 2 foreign objects, soft carry bag
100-1605U	Choker Accessory Kit: 2 objects each on a string necklace

Ambu

Ambu A/S
Baltorpbakken 13
2750 Ballerup
Denmark
T +45 72 25 20 00
F +45 72 25 20 50
ambu.com