

# KingLT(S)-D™

supraglottic airway



King LTS-D shown

## Simple. Disposable. Gastric Access.

The King LT(S)-D laryngeal tube is a disposable, simple to use alternative airway device that provides superior patient ventilation. The King LTS-D additionally allows the passage of an 18 Fr gastric tube through a separate channel.

### Simple.

- Seals in the esophagus and oropharynx to provide positive pressure ventilation.
- Single inflation port.
- Color-coded full range of sizes.
- Suction port.\*

### Disposable.

- No risk of cross-contamination.

### Gastric Access.

- Suction with up to an 18 Fr suction catheter.\*

\*King LTS-D only.

# KingLT(S)-D™

supraglottic airway

## Product Information

Product	Size 2	Size 2.5	Size 3	Size 4	Size 5
King LT-D	KLTD202	KLTD2025	KLTD203	KLTD204	KLTD205
King LTS-D	n/a	n/a	KLTS403	KLTS404	KLTS405


## Sizing Information

Size	2	2.5	3	4	5
Connector Color	Green	Orange	Yellow	Red	Purple
Patent Criteria	35-45 inches (90-115 cm) or 12-25 kg	41-51 inches (105-130 cm) or 25-35 kg	4-5 feet (122-155 cm)	5-6 feet (155-180 cm)	greater than 6 feet (>180 cm)
Cuff Pressure	60 cm H <sub>2</sub> O	60 cm H <sub>2</sub> O	60 cm H <sub>2</sub> O	60 cm H <sub>2</sub> O	60 cm H <sub>2</sub> O
KLTD O.D./I.D.	11 mm/7.5 mm	11 mm/7.5 mm	14 mm/10 mm	14 mm/10 mm	14 mm/10 mm
KLTD O.D./I.D.*	n/a	n/a	18 mm/10 mm	18 mm/10 mm	18 mm/10 mm
KLTD Cuff Volume	25-35 ml	30-40 ml	45-60 ml	60-80 ml	70-90 ml
KLTS Cuff Volume	n/a	n/a	40-55 ml	50-70 ml	60-80 ml

### King LTS-D is not available in size 2 and 2.5

**KLTD:** Maximum Size Fiberoptic Bronchoscope: 7.0 mm O.D. (size 3, 4, 5) and 4.7 (size 2, 2.5); Maximum Size Tube Exchange Catheter: 19 Fr (size 3, 4, 5) and 14 Fr (size 2, 2.5); Minimum Mouth Opening: 16 mm (size 3, 4, 5) and 12 mm (size 2, 2.5)

**KLTS:** \*Ventilation Lumen is not round, but is equivalent to a 10 mm I.D. tube; Maximum Size Tube Exchange Catheter: 19 Fr; Maximum Size Fiberoptic Bronchoscope: 6 mm OD; Minimum Mouth Opening: 20 mm



King LT-D supraglottic airway

## Kingsystems

15011 Herriman Boulevard, Noblesville, IN 46060  
 p 317.776.6823 p 800.642.5464 f 317.776.6827  
[www.Kingsystems.com](http://www.Kingsystems.com)